



# Oxygen Concentrators



[DirectSupply.net](http://DirectSupply.net)

# SAVE BIG WHEN YOU OWN YOUR OXYGEN CONCENTRATORS

## Get the Industry's Top Performers

Direct Supply offers the industry's widest selection of oxygen equipment from brands you trust. Our product experts will help you choose the best units for your staff and residents.

Contact your account manager at **1-800-634-7328** today to learn more about Direct Supply's lineup of quality respiratory equipment.

## Choose from a vast selection to meet your needs

Oxygen concentrators provide respiratory users with increased mobility and flexibility for their active lifestyles. We carry a wide selection of the top models from portable to stationary to fit your community's requirements.

### 5L Oxygen Concentrators




### 10L Oxygen Concentrators




### Portable Oxygen Concentrators




You could save your community thousands of dollars every year when you purchase oxygen concentrators rather than opting for costly rental contracts. Many communities find the purchase pays for itself in the first year!

 Calculate your potential savings today at [DirectSupply.com/Oxygen](https://DirectSupply.com/Oxygen).

You'll not only see your savings analysis, but videos on preventive maintenance and a huge selection of products to suit your residents' needs, including options from Invacare – an industry leader in respiratory products.



## 5-Year Cost Comparison: Renting vs. Owning



→ **EXAMPLE ONLY** – your savings may vary. Calculate your community's potential savings using our calculator.



Find more respiratory products on [DirectSupply.net](https://DirectSupply.net)

# 5L Oxygen Concentrators



5L concentrators are a mobile and affordable way to provide up to 5 liters of continuous oxygen per minute to residents. Key selection criteria (in addition to price point) include sound level, required preventive maintenance steps and power consumption. Most brands offer optional built-in alarms for oxygen purity and flow rate, which can help provide peace of mind when owning concentrators. The comparison chart below provides a high-level overview of the models in our offering. Proceeding pages will provide additional detail on each model, including visual guides for preventive maintenance steps.



Attendant® 5L  
Oxygen Concentrator

**PAGE 6**  
#8X899



Caire® Companion 5™

**PAGE 9**  
#C5319



DeVilbiss® Compact

**PAGE 11**  
#17121

Manufacturer	Direct Supply	Caire	DeVilbiss
<b>Unique Features</b>	<ul style="list-style-type: none"> <li>Comes standard with oxygen purity and low flow rate indicators</li> <li>Compact footprint and lightweight design allows for easy storage and transportation</li> <li>Two filters need attention only once a year under normal conditions</li> <li>Easy-to-read digital hour meter display</li> </ul>	<ul style="list-style-type: none"> <li>OPI oxygen monitor keeps track of purity levels; also available without monitor</li> <li>Only two filters to change every 2 years</li> <li>Readout display provides an error code when alarms are on</li> <li>Flow meter ball drops to visually indicate loss of pressure for resident safety</li> </ul>	<ul style="list-style-type: none"> <li>Purity sensor attached to the concentrator checks the purity level of the oxygen, allowing caregivers to feel secure that the oxygen flowing is safe for the resident</li> <li>Pressure-compensated flow meter displays accurate flow</li> <li>Turn-down technology automatically activates at 2½ lpm or less to save energy and minimize stress on compressor</li> </ul>
<b>Sound Level</b>	≤50 dB	40 dB	48 dB
<b>Power Consumption</b>	360W (average)	250W at 2 lpm, 350W at max	310W at 5 liters per minute; 275W in Turn-Down (at or below 2½ lpm)
<b>Warranty</b>	3 years or 12,000 hours of compressor usage	Full 3 years with optional 2-year extension	Full 3 years

# Attendant® 5L Oxygen Concentrator


5-Liter Oxygen Concentrator



#8X899



### Once a Year:



Clean cabinet filter; replace intake filter

### Once every 2 years:



Verify the oxygen purity, flow rate and pressure of your concentrator by using an oxygen analyzer.

### When changing residents:



Discard oxygen tubing, cannula/mask, oxygen outlet connector, and humidifier bottle. Replace cabin filter. Clean concentrator cabinet. Inspect plugs, cords and components. Check oxygen concentration.

Specifications:	
Flow Rate:	0.5 - 5 lpm
O <sub>2</sub> Purity %:	93% ± 3%*
Outlet Pressure:	5 psi ± 6.23 psi
Power Consumption:	360W (average)
Sound Level:	≤50 db
Dimensions:	21 3/16" H x 13" W x 10 1/16" D
Weight:	36 lbs.
Warranty:	3 years or 12,000 hours of compressor usage

\*0.5 - 5 LPM after running for 30 minutes

- Trusted Attendant® brand
- Easy to maintain – just 2 filters are easily accessible and need attention annually
- Low flow rate and oxygen purity indicators standard
- Alarms for power failure, sensor failure and low flow lasting > 10 minutes



# Caire® Companion 5™

5-Liter Oxygen Concentrator



## Preventive Maintenance

### Once Every 2 Years:



Check and, if necessary, change the air intake filter. We also recommend changing between residents.



Check and, if necessary, change the HEPA filter. Direct Supply suggests changing filter between residents.



Verify the oxygen purity, flow rate and pressure of your concentrator by using an oxygen analyzer.

### When changing residents:



Discard oxygen tubing, cannula/mask, oxygen outlet connector, and humidifier bottle. Replace cabin filter. Clean concentrator cabinet. Inspect plugs, cords and components. Check oxygen concentration.

### Specifications:

Flow Rate:	0.5 - 5.0 lpm
O <sub>2</sub> Purity %:	90% +5.5 / -3%
Outlet Pressure:	4.6 psi
Power Consumption:	250W at 2 lpm, 350W at max.
Sound Level:	40 dB
Dimensions:	12½"W x 13½"D x 21½"H
Weight:	36 lbs.
Warranty:	3 years*



#C5319

CAIRE®

- Maintains higher O<sub>2</sub> purity as flow rate increases
- Minimizes “exhale” noise for improved resident experience
- Reduces power consumption by 15% when running at 2 lpm or below
- Optional alarm screen displays error codes to specify technical issues
- Only two filters to change every 2 years – no weekly preventive maintenance

# DeVilbiss® Compact

5-Liter Oxygen Concentrator



## Preventive Maintenance

### Once a Week:



Clean the cabinet filter.

### Once Every 2 Years:



Verify the power-loss alarm is working. Let the concentrator run for 1 minute then unplug the power cord. An alarm should sound, indicating it is working properly.



Inspect intake filter; replace as necessary.

### Once Every 3 Years:



Check and, if necessary, change the final bacteria filter. We also recommend changing between residents.



Verify the oxygen purity, flow rate and pressure of your concentrator by using an oxygen analyzer.

### When changing residents:



Discard oxygen tubing, cannula/mask, oxygen outlet connector, and humidifier bottle. Replace cabin filter. Clean concentrator cabinet. Inspect plugs, cords and components. Check oxygen concentration.



#17121

**drive** | **DeVilbiss**  
HEALTHCARE

- Purity sensor attached to the concentrator checks the purity level of the oxygen, allowing caregivers to feel secure that the oxygen flowing is safe for the resident
- Pressure-compensated flow meter displays accurate flow
- Turndown technology automatically activates at 2½ lpm or less to save energy and minimize stress on compressor

Specifications:	
Flow Rate:	1 to 5 lpm
O <sub>2</sub> Purity %:	93% ±3% at all settings
Outlet Pressure:	8.5 psi
Power Consumption:	310W at 5 liters per minute; 275W in Turn-Down (at or below 2½ lpm)
Sound Level:	48 dB
Dimensions:	13½"W x 12"D x 24½"H
Weight:	36 lbs.
Warranty:	3 years

# 10L Oxygen Concentrators

## Drive DeVilbiss® Compact 10L 10-Liter Oxygen Concentrator

10L concentrators provide up to 10 liters of oxygen per minute and are ideal for any residents with high flow needs, or those who are prescribed over 5 liters of oxygen per minute. Key selection criteria (in addition to price point) include sound level, preventive maintenance requirements, overall size, psi and presence of unique features. Higher psi levels are typically ideal for residents with longer tubing or trach connections. Some 10L models include unique features, such as dual flow meters and outputs, which allow up to 5 lpm of oxygen to be delivered to two different residents from one machine. Please review the information in the following pages to select the solution that best fits the needs of your residents.



Drive DeVilbiss® Compact 10L  
**PAGE 13**  
#4HR15



AirSep® NewLife® Intensity 10  
**PAGE 14**  
#A9972 Single Flow, #A9974 Dual Flow

Manufacturer	Invacare	AirSep
<b>Unique Features</b>	<ul style="list-style-type: none"> <li>Assembled in the USA and compatible with cylinder transfill systems</li> <li>Lightweight with a small footprint</li> <li>Exclusive DeVilbiss OSD® (Oxygen Sensing Device), standard on every unit</li> <li>Equipped with accessible patient controls, bright LEDs, protected cannula fitting and recessed humidifier nook</li> </ul>	<ul style="list-style-type: none"> <li>High outlet pressure designed to power high-acuity respiratory accessories: nebulizers, venti masks, long oxygen tubing</li> <li>Simplified tracheostomy setup with jet nebulizer attachment</li> <li>Several alarms for added security</li> </ul>
<b>Sound Level</b>	55 dB	55 dB
<b>Power Consumption</b>	639W	590W
<b>Outlet Pressure</b>	20 +/- 1.0 psi	20 psi
<b>Warranty</b>	3 years limited	Full 3 years with optional 2-year extension

### Preventive Maintenance

#### Once a Week:



Clean the cabinet filter.

#### Once Every 2 Years:



Check and, if necessary, change the inlet compressor filter. We also recommend changing between residents.



Check and, if necessary, change the bacteria filter. We also recommend changing between residents or every 17,250 hours.



Verify the power-loss alarm is working. Let the concentrator run for 1 minute then unplug the power cord. An alarm should sound, indicating it is working properly.

#### Once Every 3 Years:



Verify the oxygen purity, flow rate and pressure of your concentrator by using an oxygen analyzer.

#### When changing residents:



Discard oxygen tubing, cannula/mask, oxygen outlet connector, and humidifier bottle. Replace cabin filter. Clean concentrator cabinet. Inspect plugs, cords and components. Check oxygen concentration.



#4HR15



- Several audible and visual alarms: high pressure, high/low flow, low oxygen, power failure
- Oxygen Sensing Device (OSD) standard
- 3-year limited warranty

Specifications:	
Flow Rate:	2 - 10 lpm
O <sub>2</sub> Purity %:	87% - 96%
O <sub>2</sub> Pressure:	20 +/- 1.0 psi
Power Consumption:	639W
Sound Level:	55 dB
Dimensions:	13½"W x 12"D x 24½"H
Weight:	42 lbs.
Warranty:	3-year limited

# AirSep® NewLife® Intensity 10

10-Liter Oxygen Concentrator




#A9974



## Preventive Maintenance

### Once a Week:



Clean and, if necessary, change the intake (gross particle) filter. We also recommend changing between residents.

### Once Every 5,000 Hours:



Check and, if necessary, change the felt filter. We also recommend changing between residents.

### Once a Year:



Verify the power-loss alarm is working. Let the concentrator run for 1 minute then unplug the power cord. An alarm should sound, indicating it is working properly.

### Once Every 2 Years:



Verify the oxygen purity, flow rate and pressure of your concentrator by using an oxygen analyzer.

### When changing residents:



Discard oxygen tubing, cannula/mask, oxygen outlet connector, and humidifier bottle. Replace cabin filter. Clean concentrator cabinet. Inspect plugs, cords and components. Check oxygen concentration.

Specifications:	
Flow Rate:	2 - 10 lpm
O <sub>2</sub> Purity %:	92% at 2 - 9 lpm, 90% at 10 lpm
O <sub>2</sub> Pressure:	20 psi
Power Consumption:	590W
Sound Level:	40 dB
Dimensions:	16½"W x 14½"D x 27½"H
Weight:	58 lbs.
Warranty:	3 years

- Available with single or dual flow meters and outputs
- Optional dual flow meter on 10 lpm models allows for two residents to receive up to 5 lpm at once
- Optional second output valve can be used for nebulizer treatments
- 20 psi output ideal for high flow needs

# Portable Oxygen Concentrators

Portable oxygen concentrators are the most lightweight and mobile O<sub>2</sub> delivery solutions available. Ranging from 4.4 to 20 lbs., these compact devices allow for AC, DC or battery power to allow oxygen therapy to be delivered while on the go. Resident activities, outside physician visits, general travel and power outages are popular reasons for using portable concentrators. Key selection criteria include flow rate (be sure to note differences between continuous and intermittent delivery options), battery life, size and weight. Please review the information on the following pages to select the solution that best fits the needs of your residents.



	Caire® Eclipse 5® PAGE 17 #2RYJY, 2RYK2	Caire® Freestyle® Comfort® PAGE 18 #6YR11, #8VK42, #8VK43	Drive DeVilbiss® iGo2 PAGE 19 #347YR, #347YT
<b>Manufacturer</b>	Caire	Caire	Drive
<b>Flow Rate</b>	0.5 to 3.0 lpm continuous, 1 - 9 pulse settings	1 - 5 pulse settings	1-5 pulse settings
<b>O<sub>2</sub> Purity %</b>	90% +5.5 / -3%	90% +5.5 / -3%	90% +4/-3%
<b>Sound Level</b>	48 dBA at 3.0 lpm continuous flow; 40 dBA at 3.0 pulse dose setting	40 dBA at pulse setting 2	37 dBA at pulse setting 2
<b>Battery Life</b>	2 hours at 2.0 lpm continuous flow	Up to 4 hours with 8-cell battery pack. Up to 8 hours with 16-cell battery pack	3.5 hrs on pulse setting 2
<b>Dimensions</b>	12¾"W x 7¼"D x 19¾"H	7¾"W x 3¼"D x 10"H	3.5"W x 8.6"D x 8.4"H
<b>Weight</b>	18.4 lbs. with battery installed	5 lbs.	4.95 lbs.
<b>Warranty</b>	3 years on concentrator, 1 year on battery and accessories	3 years on concentrator, 1 year on battery and accessories	5 years device, 3 years compressor, 1 year sieve beds and battery



# Caire® Eclipse 5®

Portable Oxygen Concentrator



## Preventive Maintenance

### Daily:



Clean and, if necessary, change the air inlet filter or once per week. Replace annually.

### Once a Year:



Replace the HEPA filter. We also recommend changing between residents.



Replace the compressor intake filter. We also recommend changing between residents.



Test the oxygen purity every time the device casing is opened, which is once annually as part of yearly preventive maintenance.

### When changing residents:



Discard oxygen tubing, cannula/mask, oxygen outlet connector, and humidifier bottle. Replace cabin filter. Clean concentrator cabinet. Inspect plugs, cords and components. Check oxygen concentration.



#2RYJY

CAIRE®

- Delivers the highest pulse dose oxygen output of any portable oxygen concentrator on the market
- Features autoSAT technology, which adjusts to the resident's respiratory rate as it increases and decreases
- Continuous flow settings 0.5 to 3.0 lpm (0.5 liter increments)

Specifications:	
Flow Rate:	0.5 - 3.0 lpm
Battery Life:	2 hours at 2 lpm continuous flow 5.1 hours at a 2 setting
O <sub>2</sub> Purity %:	90% +5.5 / -3%
O <sub>2</sub> Pressure:	5 psi
Power Consumption:	52W at 1.0 lpm continuous flow; 145W at 3.0 lpm; 45W at 1.0 pulse dose; 95W at 6.0 pulse dose
Sound Level:	45 dBA at 3.0 lpm continuous flow; 40 dBA at 3.0 pulse dose setting
Dimensions:	7 <sup>1</sup> / <sub>10</sub> "D x 12 <sup>7</sup> / <sub>10</sub> "W x 19 <sup>3</sup> / <sub>10</sub> "H
Weight:	15 lbs.
Warranty:	3 years

# Caire® Freestyle® Comfort®

Portable Oxygen Concentrator

# Drive DeVilbiss® iGo2

Portable Oxygen Concentrator



#6YR11

CAIRE®

## Preventive Maintenance

### Once a Week:



Clean 2 air intake filters. Your oxygen provider may advise cleaning more frequently depending on operating conditions.

### Every 3-6 months:



Replace the nasal cannula and tubing (every 3 months if daily use), or as recommended by the cannula manufacturer or your oxygen provider.

### When changing residents:



Discard oxygen tubing, cannula/mask, oxygen outlet connector, and humidifier bottle. Replace cabin filter. Clean concentrator cabinet. Inspect plugs, cords and components. Check oxygen concentration.

Specifications:	
Flow Rate:	1 - 5 pulse settings
O <sub>2</sub> Purity %:	90% +5.5/-3%
Sound Level:	40 dBA at pulse setting 2
Battery Life:	Up to 4 hours with 8-cell battery pack. Up to 8 hours with 16-cell battery pack
Dimensions:	7 <sup>3</sup> / <sub>16</sub> "W x 3 <sup>1</sup> / <sub>16</sub> "D x 10"H
Weight:	5 lbs.
Warranty:	3 years on concentrator, 1 year on battery and accessories

- UltraSense sensitive breath technology
- Ergonomic shape
- Color display and glow in the dark keypad
- Long battery life up to 8 hours
- Internal filters designed to last the life of the device--no replacement needed

## Preventive Maintenance

### Periodic intervals:



Clean and replace the cannula and tubing as instructed by the manufacturer and your equipment provider.

### When changing residents:



Clean and disinfect the exterior covers of the device and rechargeable battery between patient use and as needed.



#347YR,  
#347YT

drive | DeVilbiss HEALTHCARE

Specifications:	
Flow Rate:	1 - 5 pulse settings
O <sub>2</sub> Purity %:	90% +4/-3%
Sound Level:	37 dBA at setting 2
Battery Life:	3.5 hours on pulse setting 2
Dimensions:	3.5"W x 8.6"D x 8.4"H
Weight:	4.95 lbs.
Warranty:	5 years device, 3 years compressor, 1 year sieve beds and battery

- SmartDose Auto-Adjusting Conserver Technology included
- Settings 1-5 pulse dose
- Rugged design, including protective overmold
- Weighs less than 5 lbs
- Most sensitive conserver trigger on the POC market



1-800-634-7328 • [DirectSupply.net](http://DirectSupply.net)